
#### **Surfer Bather**

The Surfer Bather is a safe way of securely transferring a child in and out of the bath.







- · Fits into a bath by securing the base board suction cups firmly to the bath
- Elevates and lowers the child at the touch of a button
- · Adjustable backrest and seat

Easy to fit

• Adjustable pommel available

#### Safe to use

- · Holds the child securely in the seat during the transfer into and out of the bath
- Includes side rails and leg straps for additional safety
- Fully waterproof hand control (IP67)
- Battery operated no electric cables

### **Technical Specifications**

Max user weight	50 kg / 7.9 stone
Max user height	135 cm / 4'5"
Max height	48 cm / 19 "
Min height	8.5 cm / 3 ¼"
Seat depth	62 cm / 24"
Total length	112cm/44"
Width (max-min)	67.5cm/26.5" 37 cm / 14.5"
Total Product weight	8.7kg/19 lbs
Backrest angle adjustment	Flat to 70 degrees
Back/headrest length	56cm / 22"
Airflo Compressor	Mk3

# Simple to maintain

- · Can be removed and cleaned with ease
- · Rinse with clean water after each bath and leave it to dry inflated

## Improves wellbeing

- Perfect for parents and carers who struggle to lift a child into and out of the bath
- · Bath time is a fun part of the day for both child and parent
- · Lowers the child right to the bottom of the bath so they are completely covered in water and benefit from the therapeutic properties of bathing
- Remove from the bath when not in use to enable a more flexible environment for shared bathrooms

