



MANGAR Raiser Lifting Cushion

The Raiser Lifting Cushion is an assisted living aid specifically designed to help raise a person that is seated.

EASY TO FIT

- Will fit on any armchair and can easily be transferred from chair to chair
- Extremely lightweight and portable

SIMPLE TO MAINTAIN

- Anti-microbial material that is easy to clean
- The material is designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

SAFE TO USE

- Can lift up to 190 kg
- Easy to use handset raises and lowers cushion at the touch of a button
- Aids independence

IMPROVES WELL-BEING

- Suitable for a variety of different situations, i.e. to help with getting out of bed
- Can be used to help with getting out of a car
- Reduces the risk of musculoskeletal injury for the carer

TECHNICAL SPECIFICATIONS

Max user weight	190 kg
Inflated height	23 cm
Inflated width	43 cm
Inflated length	43 cm
Maximum inlet pressures:	0.35 bar (5spi)
Max/Min temperature	+50°C/-15°
Product weight (base unit)	2.8 kg
Airflo Compressor	12
Lifting performance per charge	up to 25 lifts



Features

1. **Made from strong polyurethane** coated nylon material
2. **Seams** joined by high frequency welding



Product code: MPCA200500



Scan me

For more information visit
www.mangarhealth.com
Email: enquiries@winnicare.uk
Phone: 0800 2800 485



WIN  CARE