



Mix and match

Maybe you just need a 2-post system over a bed, or maybe you need full room covering with a 4-post system.

No matter what your needs are, the Mobile Flex is your flexible solution, with modular parts which can be fitted and combined for your purpose.

Mobile Flex

Designed for people who need an immediate and/or temporary lifting solution that can be easily and simply fitted, used and removed. When combined with the innovative Luna ceiling lift, the mobile flex is a simple and functional lifting system and the ideal choice for either home care or institutional use.

EASE OF USE

The Mobile Flex is easy to assemble. No tools are required to adjust the height and width to provide reliable temporary solutions.

FLEXIBLE APPLICATION

Adjusting the height via gas spring is carried out in just a few seconds providing immediate use for the carer and client. The telescopic tracks give enough lifting height to provide solutions that are suitable for most rooms and most transfer situations.

For more information visit

www.winnicare.uk
Email: enquiries@winnicare.uk
Phone: 0800 2800 485

ML1155 Issue 1

LIFTING CAPACITY

The Mobile Flex provides a lifting capacity up to 275kg thus utilizing the full range of the Luna ceiling lift.

RELATED PRODUCTS

- LUNA IPX4
- LUNA TWIN
- LUNA SKY



Scan here to view our full Ceiling Lift brochure

SCAN ME



2-post system

Height and width adjustable Only takes 30 minutes to assemble



4-post system

Full adjustable, room covering and free-standing Only takes 60 minutes to assemble