


Camel Lifting Cushion

Patient
Lifting

The Camel Lifting Chair will both sit up and lift a fallen person.

Backrest gives additional support during a lift



Product code: MPCA060200

Easy to use

- Lifting cushion is inflated using an Airflo Compressor using a simple hand control
- Simple to use, requiring very little training
- Can be used independently or assisted
- Suitable for indoor and outdoor use

Simple to maintain

- Can be cleaned using a proprietary non-abrasive mild liquid cleaner.
- Easy to deflate, fold away and store – carry bag provided
- Does not need a LOLER certificate
- The new Riverseal® range of materials are designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

Safe lifting

- Lifts the frailest person up to 50 stone / 320 kg from the floor
- Backrest gives additional support and is particularly suitable for plus size patients
- Stable and supportive

Improves well-being

- Provides a safe, dignified lift for a fallen person
- Individuals with a cognitive impairment are less distressed when lifted with a Camel
- Reduces the risk of musculoskeletal injury for the carer
- Backrest gives additional support and is particularly suitable for plus size or bariatric patients

Technical Specifications

Max user weight	320 kg / 50 stone
Max height	56 cm / 22"
Min height	Flat
Seat depth	47 cm / 18 1/2"
Total length	112 cm / 44"
Width	70 cm / 27 1/2"
Total product weight	6.5 kg / 14 lbs
Back/headrest length	81 cm / 32"
Airflo Compressor	24
Lifting performance per charge	4 lifts

Mangar Health

Presteigne, Powys, Wales, LD8 2UF

t 0800 2800 485 **f** 01544 260287 **e** info@mangarhealth.com

w mangarhealth.com

ML1019 Issue 1