


ALOVA® RANGE



VMA64 model



VMA63 model



VMA62 model



KEY FEATURES

- ▶ Single-piece mattress with different hardnesses made up of:
 - 6 cm molded visco elastic foam mattress top layer: density $\geq 75 \text{ Kg/m}^3$
 - 8 cm polyether foam base: density 18 Kg/m^3 (VMA62), 28 Kg/m^3 (VMA63) and $\geq 34 \text{ Kg/m}^3$ (VMA64)
- ▶ 3 different hardnesses: medium (1.7 kPa), firm (1.9 kPa) and soft (1.5 kPa)
- ▶ Integrated heel relief
- ▶ Patient weight: 20 to 120 Kg (VMA62), 20 to 150 Kg (VMA63) and 20 to 180 Kg (VMA64)
- ▶ Dimensions: 195 x 88 x 14 cm (76.7 x 34.6 x 5.5 inches)
- ▶ Mattress weight: 14 Kg (VMA62), 15 Kg (VMA63) and 16 Kg (VMA64)
- ▶ Covers available: Promust PU, Promust PU HD and Promust CIC



KEY BENEFITS

- ▶ Visco elastic foam provides excellent pressure re-distribution properties: homogeneous immersion of the body into the foam, thus increasing the bearing surface and reducing pressure applied in areas at risk of pressure ulcers.
- ▶ Different hardnesses adapted to each body area.
- ▶ Integrated heel relief to enhance blood circulation.
- ▶ Visco elastic foam provides patient comfort.

INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT MEDIUM TO HIGH RISK OF PRESSURE ULCER DEVELOPMENT, BEDRIDDEN MORE THAN 15 HOURS A DAY.

HELP IN TREATING EXISTING PRESSURE ULCERS FROM STADE 1 TO 2 COMBINED WITH TECHNICAL AID DEVICES.

PATIENT WEIGHT: 20 TO 180 KG/ 3.1 TO 28.3 ST

- ▶ Single-piece mattress made up of molded visco elastic foam top layer on polyether foam base
- ▶ Integrated heel relief
- ▶ Removable and washable integral cover
- ▶ 3 models for 3 environments: VMA62 for homecare, VMA63 for nursing homes, VMA64 for hospitals

DESCRIPTION	ITEM REFERENCES ACCORDING TO PACKAGING			FIRE STANDARDS
	CARDBOX	PALLET BOX	ROLLED AND COMPRESSED (RC)	
ALOVA® Dermalon draw sheet	—	VMA62DH-BOX	—	EN 597-1 and 2 GPEM D1 90 classe D
ALOVA® Promust PU cover	VMA62C-HIP	VMA62-HIP (10 units)	VMA62RC1-HIP (1 unit) VMA62RC-HIP (5 units)	
	VMA62C-HIP3	VMA62-HIP3 (10 units)	—	
ALOVA® Promust PU HD cover	VMA63-HIP	VMA63-HIP-BOX (10 units)	VMA63RC1-HIP (1 unit) VMA63RC-HIP (5 units) VMA63RC1-HIP3 (1 unit) VMA63RC-HIP3 (5 units)	
	VMA63HIPH	VMA63HIPH-BOX (10 units)	VMA63RC1-HIPH3 (1 unit) VMA63RC-HIPH3 (5 units)	
ALOVA® Promust CIC cover	VMA64HIPH	VMA64HIPH-BOX (10 units)	VMA64RC1-HIPH (1 unit) VMA64RC-HIPH (5 units)	
	VMA63CIC	VMA63CIC-BOX (10 units)	—	
	VMA64CIC	VMA64CIC-BOX (10 units)	—	

Other versions available: VMA63RC-HIPH and VMA64RC-HIPH (5 units)

3 years warranty
2 years warranty in nursing homes or hospitals

Cleaning instructions:
Dermalon FR* draw sheet

Cleaning instructions:
Promust PU (HD) and Promust CIC covers



Photos not contractually binding

FP000035/01.2018

