


MEDICALYS®



2nd generation MEDICALYS bed with WINNEA boards and blanket holder



Wired remote control



Caregiver's wired control console ⁽³⁾



Connector handle with reel ⁽³⁾



WINNEA® CIC board ⁽²⁾ with wall bumper



Blanket holder ⁽³⁾



Pedal: centralised braking (Directional wheel, freewheels)



⁽¹⁾ CPR - Cardio Pulmonary Resuscitation
⁽²⁾ CIC - Cross Infection Control
Core treatment by silver ions

⁽³⁾ OPTIONS

- ▶ Metal castors Ø 100 mm
- ▶ Metal castors Ø 125 mm
- ▶ Design castors Ø 125mm
- ▶ Design castors Ø 150mm
- ▶ Removable compact bed base.
- ▶ Auto-regression back rest.
- ▶ Backup battery (3 cycles).
- ▶ Blanket holder.
- ▶ Mechanical Trendelenburg - 12° max
- ▶ Electric Reverse Trendelenburg. +16°max
- ▶ Caregiver's side control console.
- ▶ Caregiver's wired control console.
- ▶ Remote handcontrol on flexible arm.

THE BED'S PROVEN DESIGN MAKES IT IDEAL FOR MEDICAL DEPARTMENTS.

CONFIGURABILITY

- ▶ Available with mechanical Trendelenburg position, electric Reverse Trendelenburg position (armchair position), electric Reverse Trendelenburg position and mechanical Trendelenburg.

SAFETY

- ▶ Mechanical Trendelenburg which is irrespective of the bed position, can be activated with or without electric power in under 15 seconds.
- ▶ Compatible with application environments 1, 2, 3 & 5.

EASE OF USE

- ▶ The patient has direct access to the 100% electric function chair position with a single touch.

PERFORMANCE

- ▶ Chair position

DESCRIPTION	DIMENSIONS	ADJUSTABLE HEIGHT	BR	LEG REST	LIFTING CAPACITY	CASTORS	STANDARDS	WATERPROOF
2ND GENERATION MEDICALYS	Bed base 90 x 200 cm EHPAD: divisible HOSPITAL: non-divisible	37 to 86.5 cm if castors Ø 125 (EHPAD Pack) 39.5 to 89 cm if castors Ø 150 (HOSPITAL Pack)	70°	12° manual cranks Knee break 26° and 16° feet adjustment by manual crank	SWL ⁽¹⁾ Mattress 20 kg Accessories 45 kg Patient 165 kg (Environments 1 and 2) Patient 195 kg (Environments 3, 4, and 5)	metal Ø 150 mm	IEC 60601 IEC 60601-2-52 compatible with environments 3, 4 and 5	LINAK IP54

- 2 accessory holders (lengthwise)
- 2 IV stand holders (at the head of the bed)
- 2 lifting pole (at the head of the bed)

Handset, electric and mechanical parts: 5-year warranty*

Boards and accessories: 1-year warranty*

*Under normal conditions of use and excluding wearing parts - please see instructions of use

⁽¹⁾ Safe working load

Photos not contractually binding

FP000028/01.2018

