



Bed with MÉDIDOM II boards, Light Dragon Tree finish.



## ERGONOMIC ASSEMBLY AND DISMANTLING



Through rondo screws for assembling the two half bed base.



Widening of the cross-brace providing lateral stability.



Control box on variable height jack  
Reduction of vibratory and sound disturbances.



Boards are fastened using EASY MOVE.

## OPTIONS & ACCESSORIES

<sup>(1)</sup> CPR - Cardio Pulmonary Resuscitation

- ▶ Auto-regression electric CPR<sup>(1)</sup>
- ▶ Auto-regression back rest
- ▶ Electric leg rest
- ▶ Centralized braking
- ▶ Wall bumper
- ▶ EVOLUTION, ACTIV & MOBILITY AID SAM range.
- ▶ Infrared remote control
- ▶ Bed skirting on castors

THE X'PRIM3 BED HAS BEEN DESIGNED ACCORDING TO THE ERGOFIX PRINCIPLE. IT IS BASED ON A PATENTED ASSEMBLY AND DISMANTLING SYSTEM. THE X'PRIM3 CAN BE ASSEMBLED AND DISMANTLED BY ONE PERSON.

### ERGONOMICS

- ▶ Boards are fastened using the Easy- Move system.
- ▶ Back and leg rests can be assembled and dismantled without tools.
- ▶ Blue parts help with assembly.

### EASE OF USE

- ▶ Through rondo screws for assembling the half bed bases.
- ▶ Side rails are fastened by indexing.
- ▶ The jack and remote control cords are fastened on the cross-brace.

### SAFETY

- ▶ Back rest can be raised up to 75°.
- ▶ Control box on variable height jack to reduce vibratory and sound disturbances.

DESCRIPTION	DIMENSIONS	ADJUSTABLE HEIGHT	BACK REST	LEG REST	LIFTING CAPACITY	CASTORS	STANDARDS	MOTORIZATION
X'PRIM3 2 functions 3 functions RAL 7035 standard color	Divisible bed base 90 x 200 cm	34 to 84 cm - ground clearance: 13 cm	0° to 75°	12° manual crank Knee break 26° and 16° feet adjustment by manual crank	SWL <sup>(1)</sup> 170 kg Accessories 15 kg Mattress 20 kg max patient weight: 135 kg	∅ 100 mm	IEC 60601-2-52 compatible with environments 3, 4 and 5	LINAK IP66

<sup>(1)</sup> Safe working load

- 2 accessory holders (lengthwise)
- 2 IV stand holders (at the head of the bed)
- 2 lifting pole (at the head of the bed)

Handset, electric and mechanical parts: 5-year warranty\*  
Boards and accessories: 1-year warranty\*

\* Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding