


XXL X'PRESS 2.0[®]



IDEAL FOR HOME USE, THE XXL X'PRESS 2.0 BED CAN BE MOVED IN VERTICAL POSITION WITHOUT ANY TOOLS



New CA40 box
easy access to the connections.
no vibratory and sound disturbances.



New locking
backlit remote control.



adjustable height
from 35 cm to 81 cm



▶ LEG REST WITH THIGH SECTION



▶ CPR BACK REST



▶ INFRA-RED REMOTE CONTROL



THE NEW XXL XPRESS 2.0 BED HAS BEEN DESIGNED TO FACILITATE:

- ASSEMBLY AND DISMANTLING
- STORING
- CLEANING
- TRANSPORT

COMFORT

- ▶ Reduction of vibratory and sound disturbances with the new LINAK CA40 box.

ERGONOMICS

- ▶ Adjustable height from 35 to 81 cm..
- ▶ Locking backlit remote control.

DESCRIPTION	DIMENSIONS	ADJUSTABLE HEIGHT	BR	LEG REST	LIFTING CAPACITY	CASTORS	STANDARDS	MOTORIZATION
XXL X'PRESS 2 electric functions	120 x 200 cm	35 to 81 cm	75°	12° manual cranks	SWL ⁽¹⁾ 315 kg mattress 30 kg patient 270 kg accessories 15 kg	Ø 100 mm with brake without thread guard	IEC 60601-2-52 compatible with environments 1, 2, 3, 4 and 5	LINAK IP66
	140 x 200 cm							
XXL X'PRESS 3 electric functions	120 x 200 cm							
	140 x 200 cm							

- 2 accessory holders (lengthwise)
- 2 IV stand holders (at the head of the bed)
- 2 lifting pole (at the head of the bed)

Handset, electric and mechanical parts: 5-year warranty*
Boards and accessories: 1-year warranty*
*Under normal conditions of use and excluding wearing parts - please see instructions of use

⁽¹⁾ Safe working load

Photos not contractually binding

FP000010/01.2018

