

# MANGAR Patient Specific Slings

Mangar disposable clip and loop slings are manufactured by the original patient specific sling company and have been widely used and trusted in the UK for 15 years. Our Slings are available to purchase direct and are also available through SBS and NHSSC Frameworks.

**FEATURES:**

- CE marked and TUV approved
- 6 familiar colours and sizes
- Fabric strength guaranteed - load tested to 272kg/290kg
- Rigorous quality checks completed before dispatch
- It's easy to switch to the Mangar Sling - No additional training required as staff will be familiar with our sling sizing and colours
- Compatible with hoists used in hospitals and care environments

**SLING COLOURS AND SIZES**



Small	Red
Medium	Yellow
Large	Green
Large, Large	Purple
Extra-large	Blue
XX Large	Brown

Weight of sling: 0.6kg

Patient specific slings avoid cross contamination and can be disposed of after the patient is charged from your care.

The patient specific **clip slings** have been designed to be used with a four-point hoist clip system, with or without powered positioning. The clip slings will support a SWL of 272kg, regardless of the size of sling.

**NEW:** The new clip design improves safety, security, and product compatibility with leading manufacturers of hoists.



The patient specific **loop slings** are designed specifically to work in conjunction with a lifter/hoist which has spreader bars equipped with a hooking system. The loop slings will support a SWL of 290kg, regardless of the size of sling.



SCAN ME

**For more information visit**

[www.winncare.uk](http://www.winncare.uk)

Email: [enquiries@winncare.uk](mailto:enquiries@winncare.uk)

Phone: 0800 2800 485

WIN  CARE