

ELK Lifting Cushion

The ELK is a lifting cushion designed to lift from the floor.

Easy to use

- The lifting cushion is inflated using an Airflo Compressor with a simple hand control
- Simple to use, requiring very little training
- Suitable for use inside or out
- Can be used by one caregiver

Simple to maintain

- Clean using a proprietary non-abrasive mild liquid cleaner.
- Easy to deflate, fold away and store
- The material is designed to an exceptionally high standard, with extraordinary durability, strength and performance in the toughest environments.

Safe lifting

- Lifts the frailest person up to 980 lbs from the floor
- Exceptionally rigid and stable
- Can be used in confined spaces where traditional lifting equipment may not fit

Improves wellbeing

- Provides a safe, dignified lift for a fallen person
- Reduces the risk of back injury for the caregiver
- Ideal for caregivers of individuals with neurological disorders who fall regularly

Product code: MPCA070400

Technical Specifications

Max user weight	980 lbs
Max height	22"
Min height	Flat
Seat depth	8"
Width	22½"
Total product weight	8 lbs
Airflo Compressor	24
Lifting performance per charge	Up to 9 lifts

