



Product code: LAA3716

## Archimedes Bath Tub Lift

The Archimedes Bath Tub Lift will lower a person into a bath tub, on a fixed position seat, using a battery powered hand control.

### Easy to fit

- Simple to fit – separates into three pieces
- Fits further back in the bath tub to maximise leg room and avoid most bath tub handles and anti-slip coatings
- Unique 4 way quick release suckers make removing the Bath Tub Lift simple – reducing strain on your back
- Lowers within 3" of the bottom of the bath tub

### Simple to maintain

- The easy to clean steel base unit has been designed with a structure that maintains strength whilst reducing weight
- Divides into three so it can be removed and cleaned with ease

### Safe to use

- Flaps have an anti-pinch design, eliminating the risk of skin being trapped between the seat and the flap on transfer
- Only lowers a person into the bath tub when there is sufficient battery power to raise them back out again
- Battery operated – no electric cables
- Fully waterproof hand control (IP67)

### Improves wellbeing

- Bathing can reduce pain and promote a feeling of wellbeing
- Equipment will reduce the amount of moving and handling required from a carer

### Technical Specifications

Max user weight	364 lbs
Max height	17"
Min height	2 1/4"
Seat depth	17 1/2"
Total length	24"
Width (ex side flaps)	15 1/2"
Width of side flap	6 1/4"
Total product weight	22 lbs
Back/headrest length	26"

