



## Pillowlift

The Handy Pillowlift assists people who need help to rise from lying to sitting.

### Easy to fit

- Sits on top of the mattress
- With no fixtures or fittings, it can be used on any size bed
- One airline connects the Airflo to the Pillowlift. Press one button to inflate

### Simple to maintain

- For a routine clean, sponge clean with a detergent solution and dry thoroughly
- The covers can be machine washed
- Spare covers can be purchased
- The of material is designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

### Safe to use

- An ultra-lightweight inflating backrest raises a person from lying to sitting
- Backrest moulds around a person to provide support
- Inflation and deflation can be stopped at any time, allowing freedom to choose the most comfortable position
- Anti-slip pillow holder

### Improves well-being

- Stable, supportive and comfortable, a person can choose to be lifted on their back or side
- The Pillowlift sits on top of the mattress and causes no disturbance to a second person in a double bed
- Minimises moving and handling risks so reducing the risk of back injury for the carer
- It is portable and lightweight which means it can be rolled up for travelling

Product code: MPCAT20200

### Technical Specifications

Max user weight	294 lbs
Min height	Flat
Total length	31 1/2"
Width	24"
Total product weight	4 1/2 lbs
Back/headrest length	28 1/2"
Backrest angle adjustment	Flat to 55°
Airflo Compressor	12
Lifting performance per charge	Up to 20 lifts