



## Bathing Cushion

The Bathing Cushion is used to help someone who has difficulty getting in and out of the bath tub.

### Easy to fit

- Fitted and removed within seconds
- Some air can be left in the backrest for greater levels of comfort
- The design adapts to the size and shape of the bath tub, so no expensive bathroom alterations are required

### Simple to maintain

- Lightweight and portable, it can be removed and cleaned with ease
- Rinse with clean water after each bath and leave it to dry/deflated
- Compact and easy to store

### Safe to use

- Only lowers a person into the bath tub if there is enough battery power to raise them back out again

- Battery operated – no electric cables
- Fully waterproof hand control (IP67)
- The material is designed to an exceptionally high standard, with extraordinary durability, strength and performance in the toughest environments.

### Improves wellbeing

- Perfect for a person who has difficulty getting in and out of the bath tub
- The cushion deflates allowing the user to lie flat in the bath tub, fully covered in water so benefiting from the therapeutic properties of bathing
- The user can take the bathing cushion with them when travelling away
- Can be removed from the bath tub when not in use so enabling a more flexible environment for shared bathrooms

Product code: MPCA030500

### Technical Specifications

|                                |                                  |
|--------------------------------|----------------------------------|
| Max user weight                | 330 lbs                          |
| Max height                     | 15 1/4"                          |
| Min height                     | Flat                             |
| Seat depth                     | 17 1/2"                          |
| Width                          | 23"                              |
| Total product weight           | 4 1/2 lbs                        |
| Back/headrest length           | 13 1/4"                          |
| Airfo Compressor               | 12                               |
| Lifting performance per charge | Up to 8 lifts (4 bathing cycles) |

