

# Camel

The Camel Lifting Chair will both sit up and lift a fallen person.



Backrest gives additional support during a lift



## Technical Specifications

Max user weight	320 kg / 50 stone
Max height	56 cm / 22 "
Min height	Flat
Seat depth	47 cm / 18 ½ "
Total length	112 cm / 44"
Width	70 cm / 27 ½ "
Total Product weight	6.5 kg / 14 lbs
Back/headrest length	81 cm / 32"
Airflo Compressor	Plus

Product code: MPCA060200

## Easy to use

- Lifting cushion is inflated using an Airflo Compressor using a simple hand control
- Simple to use, requiring very little training
- Can be used independently or assisted
- Suitable for indoor and outdoor use

## Safe lifting

- Lifts the frailest person up to 50 stone / 320 kg from the floor
- Backrest give additional support and is particularly suitable for plus size or bariatric patients
- Stable and supportive

## Simple to maintain

- Can be cleaned using a proprietary non-abrasive mild liquid cleaner.
- Easy to deflate, fold away and store – carry bag provided
- Does not need a LOLER certificate

## Improves wellbeing

- Provides a safe, dignified lift for a fallen person
- Individuals with a cognitive impairment are less distressed when lifted with a Camel
- Reduces the risk of musculoskeletal injury for the carer
- Backrest gives additional support and is particularly suitable for plus size or bariatric patients

